

Benchmarking digital displays (monitors) for histological diagnoses: the nephropathology use case

Giorgio Cazzaniga (a), ¹ Francesco Mascadri, ¹ Stefano Marletta (a), ² Alessandro Caputo, ³ Gabriele Guidi, ⁴ Giovanni Gambaro, ⁵ Albino Eccher, ⁶ Angelo Paolo Dei Tos, ⁷ Fabio Pagni (a), ¹ Vincenzo L'Imperio (a) ¹

ABSTRACT

► Additional supplemental material is published online only. To view, please visit the journal online (https://doi.org/10.1136/jcp-2024-209418).

¹Department of Medicine and Surgery, Pathology, IRCCS Fondazione San Gerardo dei Tintori, University of Milano-Bicocca, Monza, Italy ²Division of Pathology, Humanitas Cancer Center, Catania, Italy ³Department of Medicine, Surgery and Dentistry "Scuola Medica Salernitana", San Giovanni di Dio e Ruggi d'Aragona University Hospital. University of Salerno, Salerno, Italv ⁴Medical Physics Unit, University Hospital of Modena, Modena, Italy ⁵Division of Nephrology, Department of Medicine,

Department of Medicine, University of Verona, Verona, Italy

⁶Department of Medical and Surgical Sciences for Children and Adults, University Hospital of Modena, University of Modena, Italy ⁷Surgical Pathology and Cytopathology Unit, Department of Medicine-DIMED, University of Padua School of Medicine, Padua, Italy

Correspondence to

Dr Vincenzo L'Imperio; vincenzo. limperio@gmail.com

Received 18 January 2024 Accepted 6 March 2024

Check for updates

© Author(s) (or their employer(s)) 2024. Re-use permitted under CC BY. Published by BMJ.

To cite: Cazzaniga G, Mascadri F, Marletta S, et al. J Clin Pathol Epub ahead of print: [please include Day Month Year]. doi:10.1136/ jcp-2024-209418 **Aim** The digital transformation of the pathology laboratory is being continuously sustained by the introduction of innovative technologies promoting whole slide image (WSI)-based primary diagnosis. Here, we proposed a real-life benchmark of a pathologydedicated medical monitor for the primary diagnosis of renal biopsies, evaluating the concordance between the 'traditional' microscope and commercial monitors using WSI from different scanners.

Methods The College of American Pathologists WSI validation guidelines were used on 60 consecutive renal biopsies from three scanners (Aperio, 3DHISTECH and Hamamatsu) using pathology-dedicated medical grade (MG), professional grade (PG) and consumer-off-the-shelf (COTS) monitors, comparing results with the microscope diagnosis after a 2-week washout period.

Results MG monitor was faster (1090 vs 1159 vs 1181 min, delta of 6–8%, p<0.01), with slightly better performances on the detection of concurrent diseases compared with COTS (κ =1 vs 0.96, 95% CI=0.87 to 1), but equal concordance to the commercial monitors on main diagnosis (κ =1). Minor discrepancies were noted on specific scores/classifications, with MG and PG monitors closer to the reference report (r=0.98, 95% CI=0.83 to 1 vs 0.98, 95% CI=0.83 to 1 vs 0.91, 95% CI=0.76 to 1, κ=0.93, 95% CI=077 to 1 vs 0.93, 95% CI=0.77 to 1 vs 0.86, 95% CI=0.64 to 1, κ =1 vs 0.50, 95% CI=0 to 1 vs 0.50, 95% CI=0 to 1, for IgA, antineutrophilic cytoplasmic antibody and lupus nephritis, respectively). Streamlined Pipeline for Amyloid detection through congo red fluorescence Digital Analysis detected amyloidosis on both monitors (4 of 30, 13% cases), allowing detection of minimal interstitial deposits with slight overestimation of the Amyloid Score (average 6 vs 7).

Conclusions The digital transformation needs careful assessment of the hardware component to support a smart and safe diagnostic process. Choosing the display for WSI is critical in the process and requires adequate planning.

INTRODUCTION

A fully digital transition in nephropathology requires an investment in the pathologists' workstation, whose most debated variable is represented by the monitor or screen settings.¹ Once the crucial choice between routine and medical device is made, the following validation process should be based according to the College of American Pathologists

WHAT IS ALREADY KNOWN ON THIS TOPIC

- ⇒ A fully digital transition in nephropathology requires an investment in the pathologists' workstation, whose most debated variable is represented by the monitor or screen settings.
- ⇒ Recently, pathology-oriented medical devices receiving the Food and Drug Administration (FDA) approval entered the commercial market, and preliminary attempts to benchmark these instruments for the primary diagnosis were published.
- ⇒ Commercial and FDA-approved pathologydedicated medical monitors are compared for primary diagnosis in nephropathology using College of American Pathologists whole slide images (WSIs) validation guidelines to obtain a real-world benchmark.

WHAT THIS STUDY ADDS

- ⇒ In nephropathology, primary diagnosis can be rendered faster on medical monitors, with equal accuracy as compared with commercial ones but slightly better performances on the detection of subtle/incipient secondary/ concurrent diseases.
- ⇒ Medical monitors may allow a more precise definition of prognostic scores/classification of glomerular diseases, closer to the traditional microscope evaluation.
- ⇒ All displays allow the employment of computational tools (eg, for amyloid detection, Streamlined Pipeline for Amyloid detection through congo red fluorescence Digital Analysis), with even better diagnostic performances as compared with the traditional microscope.

HOW THIS STUDY MIGHT AFFECT RESEARCH, PRACTICE OR POLICY

- ⇒ The introduction of digital pathology, independently from the workstation setting, is safe and can promote the application of innovative artificial intelligence algorithms.
- ⇒ This benchmark effort can help pathologists and stakeholders on the correct choice of the most suitable monitors/displays for primary diagnosis in digital pathology and nephropathology.

(CAP) guidelines.² Recently, pathology-oriented medical devices receiving the Food and Drug Administration (FDA) approval³ entered the commercial





Figure 1 Study design. A consecutive series of 60 renal biopsies scanned with three different devices (Aperio CS2/FL, MIDI II and NanoZoomer S60) were retrieved. Histochemical, IF and IHC glass slides were evaluated under the traditional microscope by an expert nephropathologist, and re-evaluated after a 2-week washout period through commercial and pathology-dedicated medical monitors. Review results and time required for the diagnosis were recorded and results analysed to obtain visualisation device benchmarking. COTS, consumer-off-the-shelf; DNAJB9, Dnaj (hsp40) homologue, subfamily b, member 9; IF, immunofluorescence; IHC, immunohistochemistry; MG, medical grade; PAS, periodic acid–Schiff; PG, professional grade; PLA2R, phospholipase A2 receptor; THSD7A, thrombospondin type 1 containing 7A domain.

market, and preliminary attempts to benchmark these instruments for the primary diagnosis were published.^{4 5} Here, we aim at performing a 'stress test' for FDA-approved pathologydedicated medical and different commercial monitors in the special field of nephropathology,⁶ evaluating the impact on the assessment of granular analytical variables that might affect the diagnostic act.⁷ To perform this analysis, the CAP whole slide images (WSIs) validation guidelines² were strictly followed to comprehensively review routine histochemistry, immunofluorescence (IF) and immunohistochemistry (IHC) stains to obtain a real-world benchmark of monitors on the complex renal biopsies use case.⁸

METHODS

Cases

The design of the study is reported in figure 1. A consecutive series of 60 renal biopsies with relative WSIs were retrieved from a multicentre dataset of fully anonymised cases during a PNRR study,^{6 7} as recommended by CAP guidelines.² To ensure that the complexity of the nephropathology routine was adequately represented in the validation process, all the histochemistry (H&E, Periodic acid–Schiff (PAS), Jones methenamine silver, Masson trichrome and Congo red), IF (IgG, IgA, IgM, C3, C1q, kappa and lambda light chains) and IHC stains (phospholipase A2 receptor (PLA2R)⁹; thrombospondin type 1 containing 7A domain (THSD7A)¹⁰ and DnaJ heat shock

protein family (Hsp40) member B9 (DNAJB9)¹¹) belonging to the retrieved cases were re-evaluated by the same nephropathologist after a washout period of 2 weeks. Four different visualisation systems were used and compared for the purpose (table 1):

- 1. Reference: glass slides under the traditional microscope (Olympus BX41, Shinjuku, Tokyo, Japan, for light microscopy and Zeiss AX10, Oberkochen, Germany, for IF)
- 2. Pathology-dedicated medical grade (MG) monitor (BARCO MDPC-8127, Courtrai, Belgium)
- 3. Professional grade (PG) commercial monitor (Philips 276E8VJSB/00, Amsterdam, the Netherlands)
- 4. Consumer-off-the-shelf (COTS) commercial monitor (HANNS-G HP205, Taipei, Taiwan)

WSIs were obtained from three scanning devices available (Aperio CS2/FL, Leica Biosystem, Nussloch, Germany; MIDI II, 3DHISTECH, Budapest, Hungary; NanoZoomer S60, Hamamatsu, Shizuoka, Japan). For the IF scanning process, exposure time was set manually and previews were assessed to obtain a final result as close as possible to the one observed by physical fluorescence microscope, as previously described.⁷ To account for possible speed/accuracy variability in WSI interpretation with different monitors, the displays were placed on the same desk under identical environmental (eg, light) conditions. Moreover, to evaluate the impact of different WSI navigation devices on the final time required for diagnosis, a comparative analysis

Table 1 Comparison of the technical features of traditional light microscope, commercial monitors and pathology-dedicated medical monitor						
Feature	Traditional microscope	COTS monitor	PG monitor	Pathology-dedicated MG monitor		
Screen technology	Optical lenses	LCD with LED backlight	IPS LCD with W-LED backlight	IPS LCD with LED backlight		
Active screen size	Viewing through eyepiece only	19.5" diagonal, 49.53 cm	27" diagonal, 68.6 cm	27" diagonal, 22.4×13.2"		
Resolution	Dependent on objective lens magnification	1.44 MP (1600×900 pixels) at 60 Hz	8 MP (3840×2160 pixels) at 60 Hz	8 MP (3840×2160 pixels) at 120 Hz		
Colour imaging	Colour limited by staining and light source	8-bit (16.7 million colours)	10-bit depth (1.07 billion colours)	10-bit depth (1.07 billion colours)		
Viewing angle	Limited to eyepiece field of view	170°	178°	178°		
Uniform luminance technology (ULT)	-	Not specified	Not specified	ULT		
Colour calibration	_	No professional calibration standards	N/A	sRGB, DICOM GSDF, native calibration		
Colour gamut	_	Not specified	NTSC 91%, sRGB 109%	NTSC 115%, sRGB 132%, DCI-P3 105%		
Ambient light presets	-	N/A	N/A	Yes, selectable reading room settings		
Luminance and contrast	-	250 cd/m ² ; 1000:1 contrast ratio	350 cd/m ² , 1000:1 contrast ratio	Calibrated luminance; 850 cd/m ² max; 1000:1 contrast ratio		
Input signals and ports	-	$1 \times D$ -Sub, $1 \times D$ VI-D (with HDCP)	1× DisplayPort 1.2, 2× HDMI 2.0, USB ports	2× DisplayPort 1.2, USB ports		
Environmental and safety specs	-	May not comply with medical device regulations	May not comply with medical device regulations	Compliant with multiple international safety standards (FDA 510k)		

COTS, consumer-off-the-shelf; DICOM, digital imaging and communications in medicine; DVI-D, Digital Visual Interface-Digital; FDA, Food and Drug Administration; GSDF, greyscale standard display function; HDCP, high-bandwidth Digital Content Protection; HDMI, High-Definition Multimedia Interface; IPS, in-plane switching; LCD, liquid-crystal display; LED, light-emitting diode; MG, medical grade; N/A, not applicable; PG, professional grade; USB, Universal Serial Bus; W-LED, white light-emitting diode.

was performed between conventional mouse versus integrated touchpad on the pathology-dedicated MG monitor.

Pathology review

Before starting the pathology review process, the different monitors were tested with the point of use quality assurance tool developed by the National Pathology Imaging Cooperative to prove the reliability of the visualisation chain and environmental conditions for primary diagnosis (online supplemental figure 1).¹² For each case under review, a comprehensive set of diagnostic parameters were extracted and documented (online supplemental table 1). These parameters included the main and secondary diagnoses, disease-specific scoring/classification systems (eg, for IgA nephropathy,¹³ antineutrophilic cytoplasmic antibody (ANCA)-associated glomerulonephritis¹⁴ and lupus nephritis¹⁵) and a detailed assessment of the main glomerular diagnostic parameters, including the total number of glomeruli, the count of globally and segmentally sclerotic glomeruli, and those exhibiting endocapillary and extracapillary hypercellularity. Additionally, the percentage of extent of interstitial fibrosis, tubular atrophy (IFTA) and arteriosclerosis on a scale of 0 (absent) to 3 (severe) was quantified. IF glomerular positivity and intensity (from 0 to 3+) were recorded, as well as the Amyloid Score.¹⁶ The review results, along with the time required to render a complete diagnosis for each case, have been comprehensively collected and organised into an Excel file (Microsoft, Redmond, USA). The obtained dataset is publicly available in the Bicocca Open Archive Research Data repository.¹⁷

Statistical analysis

Collected data underwent statistical analysis using Pandas and Scikit-learn Python libraries. For the evaluation of discrete variables, such as the final diagnosis, we used Cohen's kappa (κ) coefficient. For continuous variables, we applied the Pearson correlation coefficient (r). χ^2 test has been employed to assess the distribution of analysed stains and rendered diagnoses across different scanners. Analysis of variance test was used to calculate

the difference in evaluation times across the three monitors. To establish a threshold for clinical applicability, we adhered to the latest guidelines set forth by the CAP, suggesting a 95% CI as the cut-off for determining sufficient concordance.²

RESULTS

Cases

A total of 180 renal biopsies (60 for each scanner) have been re-evaluated, whose characteristics are reported in table 2. Each case consisted of a minimum set of four histochemical and eight IF stains, with occasional repetition for technical or diagnostic purposes. Some cases were subjected to specific histochemical (Congo red, n=30) or IHC stains, including PLA2R and THSD7A (n=21 each for membranous nephropathy), and DNAJB9 (n=9 for suspected fibrillary glomerulonephritis). Up to 21 different diagnosis groups were collected overall in the cohort and the consecutive enrolment allowed a random distribution of cases among the biopsies digitised with each scanner (p=0.761).

Performances on main/secondary diagnosis

The fastest method to review the whole batch of renal biopsies was the traditional microscope (990 min), followed by the MG, PG and COTS (1090 vs 1159 vs 1181 min, $\Delta 6$ –8%, p<0.01). The review on MG monitor was faster even when considering single-scanner batches (300, 370, 420 vs 307, 384, 468 vs 311, 390, 480 min, $\Delta 2$ –4%, $\Delta 4$ –5%, $\Delta 11$ –14% for 3DHISTECH, Hamamatsu and Aperio, respectively), unveiling an impact for the scanner used on the readability of the biopsies as well. The integrated touchpad slightly expedited the slide assessment, although no statistically significant impact was noted on the final diagnosis rendered on the pathology-dedicated MG monitor as compared with the conventional mouse (1090 vs 1098, p=0.094), suggesting only minor influence of the WSI navigation systems. Concordance metrics are reported in table 3. There was an optimal concordance on main diagnosis both with

		•	J .	
	Aperio CS2/ FL	3dhistech Midi II	Hamamatsu NanoZoomer S60	Total
No of cases (n)	60	60	60	180
Histochemistry (n)				
H&E	112	102	77	291
PAS	66	81	66	213
Trichrome	66	75	68	209
Jones	68	67	60	195
Congo red	5	12	13	30
IF (n)				
lgG	66	68	60	194
lgA	62	69	61	192
lgM	60	60	60	180
C3	62	63	60	185
C1q	61	65	60	186
Карра	63	72	60	195
Lambda	63	72	60	195
IHC (n)				
PLA2R	7	9	5	21
THSD7A	7	9	5	21
DNAJB9	2	1	6	9
Disease (n, %)				
MCD	8 (13)	8 (13)	6 (10)	22 (12)
MN	7 (12)	9 (15)	5 (8)	21 (12)
FSGS	7 (12)	6 (10)	7 (12)	20 (11)
IgA nephropathy	2 (3)	6 (10)	11 (18)	19 (11)
ANCA	9 (15)	4 (6)	5 (8)	18 (10)
DN	5 (8)	7 (12)	6 (10)	18 (10)
TIN	6 (10)	6 (10)	5 (8)	17 (9)
MGRS	7 (12)	4 (6)	5 (8)	16 (9)
Other	9 (15)	10 (16)	10 (16)	29 (16)

Table 2 Distribution of the available stains and final diagnosis within the case cohort retrieved, divided per scanner group

ANCA, antineutrophilic cytoplasmic antibody; DN, diabetic nephropathy; DNAJB9, Dnaj (hsp40) homologue, subfamily b, member 9; FSGS, focal segmental glomerulosclerosis; IF, immunofluorescence; IHC, immunohistochemistry; MCD, minimal change disease; MGRS, monoclonal gammopathy of renal significance; MN, membranous nephropathy; PAS, periodic acid–Schiff; PLA2R, phospholipase A2 receptor; THSD7A, thrombospondin type 1 containing 7A domain; TIN, tubulointerstitial nephritis;

pathology-dedicated MG and commercial monitors as compared with the reference microscope (κ =1), with slightly lower performances on secondary diagnosis for COTS (κ =0.96, 95% CI=0.87 to 1), reflecting one case of missed incipient light chain deposition disease (LCDD) kappa, concurrent with a prevalent light chain cast nephropathy, only detected on the microscope and on medical monitor. This minimal discrepancy can be partly explained by the comparison of IF among the three monitors, with almost perfect concordance for all antisera in the positivity assessment and a slight overestimation of the intensity with the MG versus PG and COTS monitors when compared with the microscope assessment (average 2.3+ vs 2.2+ vs 1.98+, respectively, from 2+).

Performances on scoring/classification systems

Discrepancies were noted when assessing comparability of specific scoring/classification systems as well, demonstrating slight superiority of the MG versus PG and COTS on the IgA nephropathy Oxford classification (r=0.98, 95% CI=0.83 to 1 vs 0.91, 95% CI=0.76 to 1), the

 Table 3
 Concordance between microscope versus medical device and microscope versus commercial device

and microscope versus commercial device							
	Microscope vs MG monitor	Microscope vs PG monitor	Microscope vs COTS monitor				
Main diagnosis (κ)	1	1	1				
Secondary diagnosis (κ)	1	1	0.96 (0.87 to 1)				
Score							
lgA MEST-C (r)	0.98 (0.83 to 1)	0.98 (0.83 to 1)	0.91 (0.76 to 1)				
ANCA vasculitis (κ)	0.93 (0.77 to 1)	0.93 (0.77 to 1)	0.86 (0.64 to 1)				
LN class (κ)	1	0.5 (0 to 1)	0.5 (0 to 1)				
DN class (κ)	1	1	1				
Amyloid type (κ)	1	1	1				
Glomeruli (r)							
Total	0.99 (0.85 to 1)	0.99 (0.84 to 1)	0.95 (0.80 to 1)				
Global sclerosis	1 (0.85 to 1)	0.99 (0.94 to 1)	0.94 (0.79 to 1)				
Seg sclerosis	0.99 (0.85 to 1)	0.99 (0.84 to 1)	0.68 (0.53 to 0.83)				
Endo hypercell	1 (0.85 to 1)	1 (0.85 to 1)	1 (0.85 to 1)				
Extra hypercell	1 (0.85 to 1)	0.99 (0.85 to 1)	0.83 (0.68 to 0.98)				
% IFTA (r)	0.98 (0.84 to 1)	0.98 (0.83 to 1)	0.9 (0.76 to 1)				
AS (0–3) (r)	0.95 (0.81 to 1)	0.90 (0.76 to 1)	0.78 (0.64 to 0.93)				
IF (0 vs any +) (r)							
lgG	1 (0.85 to 1)	0.97 (0.83 to 1)	0.98 (0.83 to 1)				
IgA	1 (0.85 to 1)	0.95 (0.75 to 1)	1 (0.85 to 1)				
lgM	1 (0.85 to 1)	0.94 (0.84 to 1)	1 (0.85 to 1)				
C3	1 (0.85 to 1)	0.94 (0.79 to 1)	0.85 (0.70 to 0.99)				
C1q	1 (0.85 to 1)	0.97 (0.82 to 1)	0.92 (0.78 to 1)				
Карра	0.99 (0.84 to 1)	0.96 (0.81 to 1)	0.96 (0.81 to 1)				
Lambda	1 (0.85 to 1)	0.98 (0.83 to 1)	0.96 (0.81 to 1)				

Cohen's kappa (κ) was used for discrete variables, Pearson correlation coefficient (r) for continuous ones.

AS, arteriosclerosis; COTS, consumer-off-the-shelf; DN, diabetic nephropathy; Endo hypercell, endocapillary hypercellularity; Extra hypercell, extracapillary hypercellularity; IF, immunofluorescence; IFTA, interstitial fibrosis and tubular atrophy; LN, lupus nephritis; MG, medical grade; PG, professional grade; Seg sclerosis, segmental sclerosis.

ANCA classification (κ=0.93, 95% CI=0.77 to 1 vs 0.93, 95% CI=0.77 to 1 vs 0.86, 95% CI=0.64 to 1) and the lupus nephritis classification (κ =1, vs 0.50, 95% CI=0 to 1 vs 0.50, 95% CI=0 to 1), the latter evaluated on the small batch available (six cases). This can be explained by the subanalysis on single glomerular features (total glomerular count, global and segmental sclerosis, endo/extracapillary hypercellularity), where the medical monitor outperformed the commercial ones. Notably, the COTS tended to miss an average of one glomerulus per case (average $\Delta 1.13$), along with a slight loss in the number of globally and segmentally sclerotic glomeruli, endocapillary and extracapillary hypercellularity per case (average Δ of 0.39, 0.26, 0.02, 0.12, figure 2), while the MG and PG monitors achieved comparable results with the microscope (average $\Delta < 0.1$). On the contrary, features that required a broader quantitative assessment, such as IFTA and arteriosclerosis, were not significantly affected by the switch between the medical and commercial monitors (average $\Delta < 1\%$ and average $\Delta < 0.1$).

Congo red evaluation

A perfect concordance was recorded for the interpretation of Congo red positivity with both commercial and medical monitors as compared with the traditional microscope, with 4 out of 30 cases (13%) receiving a diagnosis of amyloidosis. The application of Congo red fluorescence (CRF) as a digital alternative to the birefringence under polarised light using the Streamlined Pipeline for Amyloid detection through congo red fluorescence Digital Analysis (SPADA) pipeline¹⁶ allowed the quantification and extension assessment of amyloid deposits in the different compartments using the Amyloid Score with both monitors, with slightly higher average Amyloid Score as compared with the birefringence assessment (6 vs 7), detecting minimal interstitial deposits that were missed by polarised light in two cases.

DISCUSSION

The technological landscape is evolving rapidly, having a progressive permeation of our pathology laboratories and twisting the traditional concept of pathologists' workstation.^{18–22} In this setting, there is still a substantial lack of evidence on the potential impact of different types of displays on pathologists' diagnostic performances.²³ Currently, COTS displays are the most widely used, typically as part of an enterprise's standard workstation configuration, selected by default for their office computing



Figure 2 Comparison of the visual quality between pathologydedicated MG monitor (top, green), PG commercial monitor (middle, orange) and COTS commercial monitor (bottom, red). While the image is brighter and sharper on the MG versus PG, due to the higher luminance (850 vs 350 cd/m²), the larger screen size (27" vs 19.5") enables a superior 'panoramic' view of the slide at the same magnification (×10), accommodating more glomeruli (7 vs 5) on the MG and PG versus COTS monitors. Finally, the higher frequency rate (120 Hz vs 60 Hz) facilitates quicker evaluations and sharper images, while superior resolution and colour gamut enhance the assessment of glomerular details. COTS, consumer-off-the-shelf; MG, medical grade; PG, professional grade.

setup or included with their laptop.²⁴ However, the first MG monitors are being approved by FDA for WSI-based primary diagnosis in pathology,²⁵ and adequate benchmark efforts as well as rigorous validation tests are required to define non-inferiority of these visualisation systems over the traditional microscope diagnosis.²⁶ In the present experience, the application of CAP validation confirmed the reliability of COTS, PG and MG monitors for the primary diagnosis of renal diseases, achieving a concordance rate of >95% independently from the display or scanner used. Although the traditional microscope still remains the fastest instrument for the routine diagnosis of renal biopsies, dissecting the digital diagnostic process with the different instruments benchmarked in this study, pathology-dedicated MG monitor demonstrated highest performances on the review time required as compared with commercial/standard solutions. This can at least be partly explained by the higher refresh rate of this device as compared with the conventional commercially available ones (120 vs 60/70 Hz), thanks to the employment of two display ports and a dedicated graphic card. Moreover, the availability of additional integrated navigation devices (eg, customisable touchpad) as suitable alternatives to the conventional mouse, even if not significantly impacting on the final time required for the diagnosis, can further contribute to the simplification of WSI consultation thanks to scrolling options and pinch-to-zoom actions, as previously described.²⁷ Basically, the navigation slowness persists as a strong source of reluctance by the pathologists to switch to digital pathology, so improvements in refresh may be welcome to favour the change. Similarly, improvements in terms of quality of the images acquired with progressively new scanners impacted on the overall time required for the diagnosis, suggesting that introducing technological innovations within the pathologists' workstation can potentially reduce the hands-on time up to 14%, with repercussions on turnaround times (TATs). Slight superiority of the pathologydedicated MG monitor was noted in the detection of subtle, concurrent secondary renal diseases, as in the case of undetected incipient LCDD kappa missed on the COTS monitor, due to a more tenuous IF intensity on the commercial displays, reflecting an overall slight underestimation of the average IF intensity, as shown by the comparative results of this study. Reasons for this can be found in the highest luminance (850 vs 350 vs 250 cd/m^2 , MG vs PG vs COTS, respectively), wider colour bit depth (10 vs 8, MG/PG vs COTS, respectively) and total range of colours (1.07 billion vs 16.7 million, MG/PG vs COTS, respectively). Similarly, the availability of a wide colour gamut (sRGB 132%) enhancing the pink/violet nuances of histological preparation in the MG monitor demonstrated its role in the interpretation of the complex renal pathology histochemical stains (from H&E to the more specific PAS, Jones, trichrome and Congo red), as highlighted by the most concordant score/classification and single histological feature assessments. Moreover, the bigger size of the pathology-dedicated MG and PG monitor (27" vs 19.5") allowed a better 'panoramic' visualisation of the WSI even at lower magnifications, potentially impacting the numerical assessment of the renal structures (eg, total number of glomeruli and globally sclerosed ones), accounting for the slight differences noted in the present comparison. The switch from microscope to digital pathology and WSI enabled the application of image analysis and computational tools, as demonstrated by the amyloid use case, where the conventional assessment of birefringence under polarised light can potentially hamper the digitisation of the Congo red stain. However, the application of CRF associated with SPADA computational pipeline¹⁶ allows the automatic detection and quantification of amyloid deposits, detecting even

Original research

minimal disease involvement and showing optimal performances on both monitors. These advancements can eventually be used for automatically annotating structures for feeding/training artificial intelligence (AI) algorithms, which is potentially further facilitated by the introduction of larger, brighter screens with integrated customisable touchpads, playing a role in speeding the annotation process.

The present study demonstrated overall comparability of the benchmarked monitors on primary diagnosis, with slightly better performances of the MG/PG in ultra-specialised settings (eg, scoring systems and classifications in nephropathology). The selection of the most appropriate display for digital pathology should take into account these benefits, along with possible reduction of TATs and applicability of AI tools in the routine diagnostics, balancing the initial investments required for their purchase, which can reach up to 10 times the cost of standard consumer displays (COTS).^{5 28}

CONCLUSIONS

The digital transformation needs careful assessment of the hardware component to support a smart and safe diagnostic process. Choosing the display for WSI is critical in the process and requires adequate planning.

Handling editor Runjan Chetty.

X Alessandro Caputo @ACaputoMD and Vincenzo L'Imperio @VLimperioMD

Contributors GC and VL'I defined the study design. GC and FM performed the scanning of renal biopsies, managed the virtual slides and processed the obtained data. VL'I performed the evaluation of the biopsies. GC, SM, AC and GGu collected and analysed the data. APDT, GGa and AE provided counselling as experts in the validation of WSI and digital pathology platforms. FP and AE provided nephropathology counselling and supervision on data analysis and interpretation, as well as providing the medical and non-medical displays. VL'I, FP, GGa and AE provided the funding acquisition and administrative support. VL'I is the author responsible for the overall content as the guarantor. All authors were involved in writing the paper and had final approval of the submitted and published versions.

Funding The work has been funded by the European Union-Next Generation EU-NRRP M6C2-Investment 2.1 Enhancement and strengthening of biomedical research in the NHS (DIPLOMAT-PNRR-MR1-2022-12375735). FP includes the project in the Italian Ministry of the University MUR Dipartimenti di Eccellenza 2023-2027 (I. 232/2016, art. 1, commi 314-337). VL'I received a Ministry of Health grant (RF GR-2021-12374235 PROFOUND).

Competing interests None declared.

Patient consent for publication Not applicable.

Ethics approval Approval was obtained from the Ethics Committee of the University of Verona (PNRR-MR1-2022-12375735, 03/16/23).

Provenance and peer review Not commissioned; externally peer reviewed.

Data availability statement Data are available in a public, open access repository.

Supplemental material This content has been supplied by the author(s). It has not been vetted by BMJ Publishing Group Limited (BMJ) and may not have been peer-reviewed. Any opinions or recommendations discussed are solely those of the author(s) and are not endorsed by BMJ. BMJ disclaims all liability and responsibility arising from any reliance placed on the content. Where the content includes any translated material, BMJ does not warrant the accuracy and reliability of the translations (including but not limited to local regulations, clinical guidelines, terminology, drug names and drug dosages), and is not responsible for any error and/or omissions arising from translation and adaptation or otherwise.

Open access This is an open access article distributed in accordance with the Creative Commons Attribution 4.0 Unported (CC BY 4.0) license, which permits others to copy, redistribute, remix, transform and build upon this work for any purpose, provided the original work is properly cited, a link to the licence is given, and indication of whether changes were made. See: https://creativecommons.org/ licenses/by/4.0/.

ORCID iDs

Giorgio Cazzaniga http://orcid.org/0000-0001-5373-0611 Stefano Marletta http://orcid.org/0000-0001-7881-8767 Fabio Pagni http://orcid.org/0000-0001-9625-1637 Vincenzo L'Imperio http://orcid.org/0000-0002-9284-2998

REFERENCES

- 1 Fraggetta F, L'Imperio V, Ameisen D, et al. Best practice recommendations for the implementation of a digital pathology workflow in the anatomic pathology laboratory by the European society of digital and integrative pathology (ESDIP). *Diagnostics* (*Basel*) 2021;11:2167.
- 2 Evans AJ, Brown RW, Bui MM, et al. Validating whole slide imaging systems for diagnostic purposes in pathology. Arch Pathol Lab Med 2022;146:440–50.
- 3 Barco launches medical-grade pathology display. 2021. Available: https://www. labpulse.com/diagnostic-technologies/pathology-and-ai/pathology-histology/ article/15299674/barco-launches-medical-grade-pathology-display [Accessed 02 Dec 2023].
- 4 Abel JT, Ouillette P, Williams CL, et al. Display characteristics and their impact on digital pathology: a current review of pathologists' future 'microscope J Pathol Inform 2020;11:23.
- 5 Clarke EL, Munnings C, Williams B, et al. Display evaluation for primary diagnosis using digital pathology. J Med Imaging (Bellingham) 2020;7:027501.
- 6 L'Imperio V, Brambilla V, Cazzaniga G, et al. Digital pathology for the routine diagnosis of renal diseases: a standard model. J Nephrol 2021;34:681–8.
- 7 L'Imperio V, Casati G, Cazzaniga G, *et al.* Improvements in digital pathology equipment for renal biopsies: updating the standard model. *J Nephrol* 2023;37:221–9.
- 8 Cazzaniga G, Rossi M, Eccher A, et al. Time for a full digital approach in nephropathology: a systematic review of current artificial intelligence applications and future directions. J Nephrol 2023;37:65–76.
- 9 L'Imperio V, Pieruzzi F, Sinico RA, et al. Routine immunohistochemical staining in membranous nephropathy: in situ detection of phospholipase A2 receptor and thrombospondin type 1 containing 7A domain. J Nephrol 2018;31:543–50.
- 10 L'Imperio V, Pieruzzi FG, Sinico RA, et al. Combined plasmatic and tissue approach to membranous nephropathy-proposal of a diagnostic algorithm including immunogold labelling: changing the paradigm of a serum-based approach. Appl Immunohistochem Mol Morphol 2020;28:376–83.
- 11 L'Imperio V, Barreca A, Vergani B, et al. Destructuring glomerular diseases with structured deposits: challenges in the precision medicine era. J Nephrol 2021;34:2151–4.
- 12 Point of use QA pathology. Available: https://www.virtualpathology.leeds.ac.uk/ research/systems/pouqa/pathology/ [Accessed 11 Feb 2024].
- 13 Trimarchi H, Barratt J, Cattran DC, et al. Oxford classification of IgA nephropathy 2016: an update from the IgA nephropathy classification working group. Kidney Int 2017;91:1014–21.
- 14 Berden AE, Ferrario F, Hagen EC, et al. Histopathologic classification of ANCAassociated glomerulonephritis. JAm Soc Nephrol 2010;21:1628–36.
- 15 Bajema IM, Wilhelmus S, Alpers CE, et al. Revision of the International society of nephrology/renal pathology society classification for lupus nephritis: clarification of definitions, and modified National Institutes of health activity and chronicity indices. *Kidney Int* 2018;93:789–96.
- 16 Cazzaniga G, Bolognesi MM, Stefania MD, et al. Congo red staining in digital pathology: the streamlined pipeline for amyloid detection through Congo red fluorescence digital analysis. *Lab Invest* 2023;103:100243.
- 17 Cazzaniga G, L'Imperio V. Benchmarking digital displays (monitors) for histological diagnoses: the nephropathology use case. *Bicocca Open Archive Research Data* 2024.
- 18 Kiran N, Sapna F, Kiran F, et al. Digital pathology: transforming diagnosis in the digital age. Cureus 2023.
- 19 Marletta S, Pantanowitz L, Malvi D, et al. Validation of portable tablets for transplant pathology diagnosis according to the college of American Pathologists guidelines. Acad Pathol 2022;9:100047.
- 20 Santonicco N, Marletta S, Pantanowitz L, et al. Impact of mobile devices on cancer diagnosis in cytology. Diagn Cytopathol 2022;50:34–45.
- 21 Brunelli M, Beccari S, Colombari R, et al. iPathology cockpit diagnostic station: validation according to college of American Pathologists Pathology and laboratory quality center recommendation at the hospital trust and University of Verona. *Diagn Pathol* 2014;9 Suppl 1:S12.
- 22 Eccher A, Fontanini G, Fusco N, *et al.* Digital slides as an effective tool for programmed death ligand 1 combined positive score assessment and training: lessons learned from the 'programmed death ligand 1 key learning program in head-and-neck squamous cell carcinoma. *J Pathol Inform* 2021;12:1.
- 23 Zarella MD, Bowman D, Aeffner F, *et al*. A practical guide to whole slide imaging: a white paper from the digital pathology association. *Arch Pathol Lab Med* 2019;143:222–34.
- 24 Norgan AP, Suman VJ, Brown CL, et al. Comparison of a medical-grade monitor vs commercial off-the-shelf display for mitotic figure enumeration and small object (Helicobacter Pylori) detection. Am J Clin Pathol 2018;149:181–5.
- 25 MDPC-8127. Available: https://www.barco.com/en/product/mdpc-8127 [Accessed 02 Dec 2023].

Original research

- 26 Randell R, Ambepitiya T, Mello-Thoms C, *et al*. Effect of display resolution on time to diagnosis with virtual pathology slides in a systematic search task. *J Digit Imaging* 2015;28:68–76.
- 27 Molin J, Lundström C, Fjeld M. A comparative study of input devices for digital slide navigation. J Pathol Inform 2015;6:7.
- 28 Picture Archiving and communication systems (PACS) and guidelines on diagnostic display devices, third edition. Available: https://www.rcr.ac.uk/our-services/all-ourpublications/clinical-radiology-publications/picture-archiving-and-communicationsystems-pacs-and-guidelines-on-diagnostic-display-devices-third-edition/ [Accessed 14 Feb 2024].